

YTHDC2 Polyclonal Antibody

Catalog No: #28800

Package Size: #28800-1 50ul #28800-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	YTHDC2 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of human YTHDC2 (NP_073739.3).
Other Names	YTHDC2;CAHL
Accession No.	Uniprot:Q9H6S0GeneID:64848
Uniprot	Q9H6S0
GeneID	64848
Calculated MW	160kDa
SDS-PAGE MW	160KDa
Formulation	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.

Application Details

WB□1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using YTHDC2 Rabbit pAb.

Background

This gene encodes a member of the DEAH (Asp-Glu-Ala-His) subfamily of proteins, part of the DEAD (Asp-Glu-Ala-Asp) box family of RNA helicases. The encoded protein binds to N6-methyladenosine, a common modified RNA nucleotide that is enriched in the stop codons and 3' UTRs of eukaryotic messenger RNAs. Binding of proteins to this modified nucleotide may regulate mRNA translation and stability. This gene may be

associated with susceptibility to pancreatic cancer in human patients, and knockdown of this gene resulted in reduced proliferation in a human liver cancer cell line. [provided by RefSeq, Sep 2016]

Note: This product is for in vitro research use only