

FOXN3 Polyclonal Antibody

Catalog No: #28808

Package Size: #28808-1 50ul #28808-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	FOXN3 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Mouse
Immunogen Description	Recombinant fusion protein of human FOXN3 (NP_005188.2).
Other Names	FOXN3; C14orf116; CHES1; PRO1635; forkhead box N3
Accession No.	Swiss-Prot#:O00409NCBI Gene ID:1112
Uniprot	O00409
GeneID	1112;
Calculated MW	54kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB □ 1:500 - 1:2000

Images



Western blot analysis of extracts of mouse pancreas, using FOXN3 at 1:1000 dilution.

Background

This gene is a member of the forkhead/winged helix transcription factor family. Checkpoints are eukaryotic DNA damage-inducible cell cycle arrests at G1 and G2. Checkpoint suppressor 1 suppresses multiple yeast checkpoint mutations including *mec1*, *rad9*, *rad53* and *dun1* by activating a MEC1-independent checkpoint pathway. Alternative splicing is observed at the locus, resulting in distinct isoforms.

Note: This product is for in vitro research use only