

NDUFC2 Polyclonal Antibody

Catalog No: #28817



Package Size: #28817-1 50ul #28817-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	NDUFC2 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IHC
Species Reactivity	Human,Mouse,Rat
Immunogen Description	A synthetic peptide of human NDUFC2 (NP_004540.1).
Other Names	NDUFC2; B14.5b; CI-B14.5b; HLC-1; NADHDH2; NADH:ubiquinone oxidoreductase subunit C2
Accession No.	Swiss-Prot#:O95298NCBI Gene ID:4718
Uniprot	O95298
GeneID	4718;
Calculated MW	14kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C


Application Details

WB 1:500 - 1:2000IHC 1:100 - 1:200

Images



Western blot analysis of extracts of various cell lines, using NDUFC2 at 1:1000 dilution.



Western blot analysis of extracts from normal (control) and NDUFC2 knockout (KO) HeLa cells, using NDUFC2 at 1:1000 dilution.

Background

Note: This product is for in vitro research use only