

HHAT Polyclonal Antibody

Catalog No: #28855



Package Size: #28855-1 50ul #28855-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	HHAT Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Mouse,Rat
Immunogen Description	A synthetic peptide of human HHAT (NP_060664.2).
Other Names	HHAT; MART2; SKI1; Skn; hedgehog acyltransferase
Accession No.	Swiss-Prot#:Q5VTY9NCBI Gene ID:55733
Uniprot	Q5VTY9
GeneID	55733;
Calculated MW	63kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using HHAT at 1:1000 dilution.

Background

Note: This product is for in vitro research use only