

## AEN Polyclonal Antibody

Catalog No: #28859



Package Size: #28859-1 50ul #28859-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                       |                                                                   |
|-----------------------|-------------------------------------------------------------------|
| Product Name          | AEN Polyclonal Antibody                                           |
| Host Species          | Rabbit                                                            |
| Clonality             | Polyclonal                                                        |
| Isotype               | IgG                                                               |
| Purification          | Affinity purification                                             |
| Applications          | WB                                                                |
| Species Reactivity    | Human,Mouse,Rat                                                   |
| Immunogen Description | Recombinant fusion protein of human AEN (NP_073604.3).            |
| Other Names           | AEN; ISG20L1; pp12744; apoptosis-enhancing nuclease               |
| Accession No.         | Swiss-Prot#:Q8WTP8NCBI Gene ID:64782                              |
| Uniprot               | Q8WTP8                                                            |
| GeneID                | 64782;                                                            |
| Calculated MW         | Refer to figures                                                  |
| Formulation           | Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4. |
| Storage               | Store at -20°C                                                    |

## Application Details

WB□1:500 - 1:2000

## Background

Note: This product is for in vitro research use only