

LRRC61 Polyclonal Antibody

Catalog No: #28861



Package Size: #28861-1 50ul #28861-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	LRRC61 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Mouse,Rat
Immunogen Description	Recombinant fusion protein of human LRRC61 (NP_076431.1).
Other Names	LRRC61; HSPC295; leucine rich repeat containing 61
Accession No.	Swiss-Prot#:Q9BV99NCBI Gene ID:65999
Uniprot	Q9BV99
GeneID	65999;
Calculated MW	Refer to figures
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB□1:500 - 1:2000

Background

Note: This product is for in vitro research use only