

FAM83A Polyclonal Antibody

Catalog No: #28869



Package Size: #28869-1 50ul #28869-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	FAM83A Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of human FAM83A (NP_116288.2).
Other Names	FAM83A; BJ-TSA-9; protein FAM83A
Accession No.	Swiss-Prot#:Q86UY5NCBI Gene ID:84985
Uniprot	Q86UY5
GeneID	84985;
Calculated MW	47kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using FAM83A at 1:1000 dilution.

Background

Note: This product is for in vitro research use only