

GIMAP8 Polyclonal Antibody

Catalog No: #28872

Package Size: #28872-1 50ul #28872-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	GIMAP8 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of human GIMAP8 (NP_783161.1).
Other Names	GIMAP8; IAN-9; IAN6; IAN9; IANT; GTPase IMAP family member 8
Accession No.	Swiss-Prot#:Q8ND71NCBI Gene ID:155038
Uniprot	Q8ND71
GeneID	155038;
Calculated MW	75kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using GIMAP8 at 1:1000 dilution.

Background

This gene encodes a protein belonging to the GTP-binding superfamily and to the immuno-associated nucleotide (IAN) subfamily of nucleotide-binding proteins. In humans, the IAN subfamily genes are located in a cluster at 7q36.1.

Note: This product is for in vitro research use only