

QSOX2 Polyclonal Antibody

Catalog No: #28873



Package Size: #28873-1 50ul #28873-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	QSOX2 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Mouse
Immunogen Description	Recombinant fusion protein of human QSOX2 (NP_859052.3).
Other Names	QSOX2; QSCN6L1; SOXN; sulfhydryl oxidase 2
Accession No.	Swiss-Prot#:Q6ZRP7NCBI Gene ID:169714
Uniprot	Q6ZRP7
GeneID	169714;
Calculated MW	85kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using QSOX2 at 1:1000 dilution.

Background

Note: This product is for in vitro research use only