

PEX2 Polyclonal Antibody

Catalog No: #28908



Package Size: #28908-1 50ul #28908-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	PEX2 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IHC
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of human PEX2 (NP_001165558.1).
Other Names	PEX2;PAF1;PBD5A;PBD5B;PMP3;PMP35;PXMP3;RNF72;ZWS3
Accession No.	Uniprot:P28328GenelD:5828
Uniprot	P28328
GenelD	5828
Calculated MW	35kDa
SDS-PAGE MW	35kDa
Formulation	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.


Application Details

WB1:200 - 1:2000IHC1:50 - 1:200

Images



Western blot analysis of extracts of various cell lines, using PEX2 antibody.



Immunohistochemistry of paraffin-embedded human placenta using PEX2 Rabbit pAb.

Background

This gene encodes an integral peroxisomal membrane protein required for peroxisome biogenesis. The protein is thought to be involved in peroxisomal matrix protein import. Mutations in this gene result in one form of Zellweger syndrome and infantile Refsum disease. Alternative splicing results in multiple transcript variants encoding the same protein.

Note: This product is for in vitro research use only