

STYX Polyclonal Antibody

Catalog No: #28913



Package Size: #28913-1 50ul #28913-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	STYX Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human
Immunogen Description	Recombinant fusion protein of human STYX (NP_660294.1).
Other Names	STYX; serine/threonine/tyrosine-interacting protein
Accession No.	Swiss-Prot#:Q8WUJ0NCBI Gene ID:6815
Uniprot	Q8WUJ0
GeneID	6815;
Calculated MW	25kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:200 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using STYX at 1:1000 dilution.

Background

The protein encoded by this gene is a pseudophosphatase, able to bind potential substrates but lacking an active catalytic loop. The encoded protein may be involved in spermiogenesis. Two transcript variants encoding the same protein have been found for these genes.

Note: This product is for in vitro research use only