

ZSCAN21 Polyclonal Antibody

Catalog No: #28917



Package Size: #28917-1 50ul #28917-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	ZSCAN21 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of human ZSCAN21 (NP_666019.1).
Other Names	ZSCAN21; NY-REN-21; ZNF38; Zipro1; zinc finger and SCAN domain containing 21
Accession No.	Swiss-Prot#:Q9Y5A6NCBI Gene ID:7589
Uniprot	Q9Y5A6
GeneID	7589;
Calculated MW	54kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:200 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using ZSCAN21 at 1:1000 dilution.

Background

Note: This product is for in vitro research use only