

MTMR6 Polyclonal Antibody

Catalog No: #28927



Package Size: #28927-1 50ul #28927-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	MTMR6 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of human MTMR6 (NP_004676.3).
Other Names	MTMR6; myotubularin-related protein 6
Accession No.	Swiss-Prot#:Q9Y217NCBI Gene ID:9107
Uniprot	Q9Y217
GeneID	9107;
Calculated MW	72kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:200 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using MTMR6 at 1:1000 dilution.

Background

Note: This product is for in vitro research use only