

GFPT2 Polyclonal Antibody

Catalog No: #28935



Package Size: #28935-1 50ul #28935-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	GFPT2 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IF
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of human GFPT2 (NP_005101.1).
Other Names	GFPT2; GFAT; GFAT 2; GFAT2; glutamine-fructose-6-phosphate transaminase 2
Accession No.	Swiss-Prot#:O94808NCBI Gene ID:9945
Uniprot	O94808
GeneID	9945;
Calculated MW	77kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C


Application Details

WB 1:200 - 1:2000IF 1:50 - 1:200

Images



Western blot analysis of extracts of various cell lines, using GFPT2 at 1:1000 dilution.



Immunofluorescence analysis of L929 cells using GFPT2 at dilution of 1:100. Blue: DAPI for nuclear staining.

Background

Note: This product is for in vitro research use only