

FNDC3A Polyclonal Antibody

Catalog No: #28952



Package Size: #28952-1 50ul #28952-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	FNDC3A Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human, Mouse
Immunogen Description	Recombinant fusion protein of human FNDC3A (NP_001265367.1).
Other Names	FNDC3A; FNDC3; HUGO; bA203I16.1; bA203I16.5; fibronectin type III domain containing 3A
Accession No.	Swiss-Prot#: Q9Y2H6 NCBI Gene ID: 22862
Uniprot	Q9Y2H6
GeneID	22862;
Calculated MW	Refer to figures
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB □ 1:200 - 1:2000

Background

Note: This product is for in vitro research use only