

## MORC3 Polyclonal Antibody

Catalog No: #28955



Package Size: #28955-1 50ul #28955-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	MORC3 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human, Mouse
Immunogen Description	Recombinant fusion protein of human MORC3 (NP_056173.1).
Other Names	MORC3; NXP2; ZCW5; ZCWCC3; MORC family CW-type zinc finger 3
Accession No.	Swiss-Prot#:Q14149NCBI Gene ID:23515
Uniprot	Q14149
GeneID	23515;
Calculated MW	Refer to figures
Formulation	Avoid freeze / thaw cycles.   Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

## Application Details

WB □ 1:200 - 1:2000

## Background

This gene encodes a protein that localizes to the nuclear matrix and forms nuclear bodies via an ATP-dependent mechanism. The protein is predicted to have coiled-coil and zinc finger domains and has RNA binding activity. Alternative splicing produces multiple transcript variants encoding distinct isoforms.

Note: This product is for in vitro research use only