

# TSPAN17 Polyclonal Antibody

Catalog No: #28958

Package Size: #28958-1 50ul #28958-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

|                       |                                                                   |
|-----------------------|-------------------------------------------------------------------|
| Product Name          | TSPAN17 Polyclonal Antibody                                       |
| Host Species          | Rabbit                                                            |
| Clonality             | Polyclonal                                                        |
| Isotype               | IgG                                                               |
| Purification          | Affinity purification                                             |
| Applications          | WB                                                                |
| Species Reactivity    | Mouse                                                             |
| Immunogen Description | Recombinant fusion protein of human TSPAN17 (NP_001006617.2).     |
| Other Names           | TSPAN17; FBX23; FBXO23; TM4SF17; tetraspanin-17                   |
| Accession No.         | Swiss-Prot#:Q96FV3NCBI Gene ID:26262                              |
| Uniprot               | Q96FV3                                                            |
| GeneID                | 26262;                                                            |
| Calculated MW         | 35kDa                                                             |
| Formulation           | Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4. |
| Storage               | Store at -20°C                                                    |

## Application Details

WB 1:200 - 1:2000

## Images



Western blot analysis of extracts of mouse stomach, using TSPAN17 at 1:1000 dilution.

## Background

This gene encodes a member of the transmembrane 4 superfamily. It is characterized by four tetraspanin transmembrane segments. The function of this gene has not yet been determined.

---

Note: This product is for in vitro research use only