

OXGR1 Polyclonal Antibody

Catalog No: #28963



Package Size: #28963-1 50ul #28963-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	OXGR1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Mouse,Rat
Immunogen Description	A synthetic peptide of human OXGR1 (NP_543008.3).
Other Names	OXGR1; GPR80; GPR99; P2RY15; P2Y15; aKGR; oxoglutarate receptor 1
Accession No.	Swiss-Prot#:Q96P68NCBI Gene ID:27199
Uniprot	Q96P68
GeneID	27199;
Calculated MW	38kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:200 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using OXGR1 at 1:1000 dilution.

Background

Note: This product is for in vitro research use only