

THEG Polyclonal Antibody

Catalog No: #28969



Package Size: #28969-1 50ul #28969-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	THEG Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of human THEG (NP_057669.1).
Other Names	THEG; CT56; THEG1; theg spermatid protein
Accession No.	Swiss-Prot#:Q9P2T0NCBI Gene ID:51298
Uniprot	Q9P2T0
GeneID	51298;
Calculated MW	Refer to figures
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB□1:200 - 1:2000

Background

This gene is specifically expressed in the nucleus of haploid male germ cells. The orthologous gene in mice encodes a protein that may play a role in protein assembly through interactions with T-complex protein 1 subunit epsilon. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene.

Note: This product is for in vitro research use only