

LRRC4C Polyclonal Antibody

Catalog No: #28997



Package Size: #28997-1 50ul #28997-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	LRRC4C Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Rat
Immunogen Description	Recombinant fusion protein of human LRRC4C (NP_065980.1).
Other Names	LRRC4C; NGL-1; NGL1; leucine rich repeat containing 4C
Accession No.	Swiss-Prot#:Q9HCJ2NCBI Gene ID:57689
Uniprot	Q9HCJ2
GeneID	57689;
Calculated MW	72kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:200 - 1:2000

Images



Western blot analysis of extracts of rat eye, using LRRC4C at 1:1000 dilution.

Background

Note: This product is for in vitro research use only