

RBM25 Polyclonal Antibody

Catalog No: #29307



Package Size: #29307-1 50ul #29307-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	RBM25 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IF
Species Reactivity	Human,Mouse
Immunogen Description	Recombinant fusion protein of human RBM25 (NP_067062.1).
Other Names	RBM25; NET52; RED120; RNPC7; S164; Snu71; fSAP94; RNA-binding protein 25
Accession No.	Swiss-Prot#:P49756NCBI Gene ID:58517
Uniprot	P49756
GeneID	58517;
Calculated MW	120kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C


Application Details

WB 1:200 - 1:2000IF 1:50 - 1:200

Images



Western blot analysis of extracts of various cell lines, using RBM25 at 1:1000 dilution.



Immunofluorescence analysis of U-2 OS cells using RBM25 at dilution of 1:100. Blue: DAPI for nuclear staining.

Background

Note: This product is for in vitro research use only