

DPEP2 Polyclonal Antibody

Catalog No: #29311



Package Size: #29311-1 50ul #29311-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	DPEP2 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of human DPEP2 (NP_071750.1).
Other Names	DPEP2; MBD2; dipeptidase 2
Accession No.	Swiss-Prot#:Q9H4A9NCBI Gene ID:64174
Uniprot	Q9H4A9
GeneID	64174;
Calculated MW	53kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:200 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using DPEP2 at 1:1000 dilution.

Background

Note: This product is for in vitro research use only