

## DCAF11 Polyclonal Antibody

Catalog No: #29318

Package Size: #29318-1 50ul #29318-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	DCAF11 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of human DCAF11 (NP_079506.3).
Other Names	DCAF11; GL014; PRO2389; WDR23; DDB1 and CUL4 associated factor 11
Accession No.	Swiss-Prot#:Q8TEB1NCBI Gene ID:80344
Uniprot	Q8TEB1
GeneID	80344;
Calculated MW	61kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

## Application Details

WB 1:200 - 1:2000

## Images



Western blot analysis of extracts of various cell lines, using DCAF11 at 1:1000 dilution.

## Background

This gene encodes a WD repeat-containing protein that interacts with the COP9 signalosome, a macromolecular complex that interacts with cullin-RING E3 ligases and regulates their activity by hydrolyzing cullin-Nedd8 conjugates. Multiple alternatively spliced transcript variants have been found for this gene.

---

Note: This product is for in vitro research use only