

UBXN6 Polyclonal Antibody

Catalog No: #29319



Package Size: #29319-1 50ul #29319-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	UBXN6 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human, Mouse
Immunogen Description	Recombinant fusion protein of human UBXN6 (NP_001164562.1).
Other Names	UBXN6; UBXD1; UBXDC2; UBX domain protein 6
Accession No.	Swiss-Prot#: Q9BZV1 NCBI Gene ID: 80700
Uniprot	Q9BZV1
GeneID	80700;
Calculated MW	Refer to figures
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB □ 1:200 - 1:2000

Background

Note: This product is for in vitro research use only