

SLIRP Polyclonal Antibody

Catalog No: #29320



Package Size: #29320-1 50ul #29320-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	SLIRP Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IHC
Species Reactivity	Human,Mouse
Immunogen Description	A synthetic peptide of human SLIRP (NP_112487.1).
Other Names	SLIRP;C14orf156;DC50;PD04872
Accession No.	Uniprot:Q9GZT3GeneID:81892
Uniprot	Q9GZT3
GeneID	81892
Calculated MW	12kDa
SDS-PAGE MW	12kDa
Formulation	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.


Application Details

WB 1:200 - 1:2000IHC 1:50 - 1:200

Images



Western blot analysis of extracts of various cell lines, using SLIRP antibody.



Immunohistochemistry of paraffin-embedded Human colon carcinoma using SLIRP Rabbit pAb.

Background

Steroid receptor RNA activator (SRA, or SRA1; MIM 603819) is a complex RNA molecule containing multiple stable stem-loop structures that functions in coactivation of nuclear receptors. SLIRP interacts with stem-loop structure-7 of SRA (STR7) and modulates nuclear receptor transactivation (Hatchell et al., 2006 [PubMed 16762838]).

Note: This product is for in vitro research use only