

PARP9 Polyclonal Antibody

Catalog No: #29322



Package Size: #29322-1 50ul #29322-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	PARP9 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human
Immunogen Description	Recombinant fusion protein of human PARP9 (NP_001139577.1).
Other Names	PARP9; ARTD9; BAL; BAL1; MGC:7868; poly [ADP-ribose] polymerase 9
Accession No.	Swiss-Prot#:Q8IXQ6NCBI Gene ID:83666
Uniprot	Q8IXQ6
GeneID	83666;
Calculated MW	Refer to figures
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB□1:200 - 1:2000

Background

Note: This product is for in vitro research use only