

TRAPPC9 Polyclonal Antibody

Catalog No: #29323

Package Size: #29323-1 50ul #29323-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	TRAPPC9 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Mouse
Immunogen Description	Recombinant fusion protein of human TRAPPC9 (NP_113654.4).
Other Names	TRAPPC9; IBP; IKBKBBP; MRT13; NIBP; T1; TRS120; trafficking protein particle complex 9
Accession No.	Swiss-Prot#:Q96Q05NCBI Gene ID:83696
Uniprot	Q96Q05
GeneID	83696;
Calculated MW	150kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:200 - 1:2000

Images



Western blot analysis of extracts of mouse pancreas, using TRAPPC9 at 1:1000 dilution.

Background

This gene encodes a protein that likely plays a role in NF-kappa-B signaling. Mutations in this gene have been associated with autosomal-recessive mental retardation. Alternatively spliced transcript variants have been described.

Note: This product is for in vitro research use only