

NICN1 Polyclonal Antibody

Catalog No: #29325

Package Size: #29325-1 50ul #29325-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	NICN1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Mouse,Rat
Immunogen Description	Recombinant fusion protein of human NICN1 (NP_115692.1).
Other Names	NICN1; nicolin-1
Accession No.	Swiss-Prot#:Q9BSH3NCBI Gene ID:84276
Uniprot	Q9BSH3
GeneID	84276;
Calculated MW	18kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:200 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using NICN1 at 1:1000 dilution.

Background

This protein encoded by this gene localizes to the nucleus and is expressed in numerous tissues including brain, testis, liver, and kidney. This refseq contains genomic sequence in its 3' UTR which is not supported by experimental evidence. Computer predictions indicate that this region of the 3' UTR contains hairpin-forming self-complementary sequence which is possibly excised after transcription. This gene has a pseudogene on chromosome X.

Note: This product is for in vitro research use only