

TRAPPC6B Polyclonal Antibody

Catalog No: #29341



Package Size: #29341-1 50ul #29341-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	TRAPPC6B Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IHC
Species Reactivity	Mouse,Rat
Immunogen Description	Recombinant fusion protein of human TRAPPC6B (NP_803235.1).
Other Names	TRAPPC6B;TPC6
Accession No.	Uniprot:Q86SZ2GeneID:122553
Uniprot	Q86SZ2
GeneID	122553
Calculated MW	18kDa
SDS-PAGE MW	18kDa
Formulation	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.


Application Details

WB 1:200 - 1:2000IHC 1:50 - 1:200


Images



Western blot analysis of extracts of various cell lines, using TRAPPC6B antibody.



Immunohistochemistry of paraffin-embedded Rat brain using TRAPPC6B Rabbit pAb.



Immunohistochemistry of paraffin-embedded Mouse brain using TRAPPC6B Rabbit pAb.

Background

TRAPPC6B is a component of TRAPP complexes, which are tethering complexes involved in vesicle transport (Kummel et al., 2005 [PubMed 16025134]).[supplied by OMIM, Mar 2008]

Note: This product is for in vitro research use only