

RNF149 Polyclonal Antibody

Catalog No: #29369



Package Size: #29369-1 50ul #29369-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	RNF149 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IHC,IF
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of human RNF149 (NP_775918.2).
Other Names	RNF149; DNATP2; ring finger protein 149
Accession No.	Swiss-Prot#:Q8NC42NCBI Gene ID:284996
Uniprot	Q8NC42
GeneID	284996;
Calculated MW	43kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C


Application Details

WB 1:500 - 1:2000IHC 1:50 - 1:200IF 1:50 - 1:200


Images




Western blot analysis of extracts of various cell lines, using RNF149 at 1:1000 dilution.




Immunofluorescence analysis of C6 cells using RNF149 Polyclonal at dilution of 1:100 (40x lens). Blue: DAPI for nuclear staining.



Immunohistochemistry of paraffin-embedded rat brain using RNF149 at dilution of 1:100 (40x lens).



Immunohistochemistry of paraffin-embedded human stomach using RNF149 at dilution of 1:100 (40x lens).



Immunohistochemistry of paraffin-embedded mouse kidney using RNF149 at dilution of 1:100 (40x lens).

Background

Note: This product is for in vitro research use only