

MAN2C1 Polyclonal Antibody

Catalog No: #29387



Package Size: #29387-1 50ul #29387-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	MAN2C1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Mouse
Immunogen Description	Recombinant fusion protein of human MAN2C1 (NP_006706.2).
Other Names	MAN2C1; MAN6A8; MANA; MANA1; alpha-mannosidase 2C1
Accession No.	Swiss-Prot#:Q9NTJ4NCBI Gene ID:4123
Uniprot	Q9NTJ4
GeneID	4123;
Calculated MW	128kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using MAN2C1 at 1:1000 dilution.

Background

Note: This product is for in vitro research use only