

TTC3 Polyclonal Antibody

Catalog No: #29401



Package Size: #29401-1 50ul #29401-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	TTC3 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of human TTC3 (NP_003307.3).
Other Names	TTC3; DCRR1; RNF105; TPRDIII; E3 ubiquitin-protein ligase TTC3
Accession No.	Swiss-Prot#:P53804NCBI Gene ID:7267
Uniprot	P53804
GeneID	7267;
Calculated MW	280kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using TTC3 at 1:1000 dilution.

Background

Note: This product is for in vitro research use only