

CLSTN3 Polyclonal Antibody

Catalog No: #29411

Package Size: #29411-1 50ul #29411-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	CLSTN3 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Mouse
Immunogen Description	Recombinant fusion protein of human CLSTN3 (NP_055533.2).
Other Names	CLSTN3; CDHR14; CSTN3; alcbeta; calsyntenin-3
Accession No.	Swiss-Prot#:Q9BQT9NCBI Gene ID:9746
Uniprot	Q9BQT9
GeneID	9746;
Calculated MW	105kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of mouse brain, using CLSTN3 at 1:1000 dilution.

Background

Note: This product is for in vitro research use only