

PLL_P Polyclonal Antibody

Catalog No: #29445



Package Size: #29445-1 50ul #29445-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	PLL _P Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IF
Species Reactivity	Mouse,Rat
Immunogen Description	A synthetic peptide of human PLL _P (NP_057077.1).
Other Names	PLL _P ; PMLP; TM4SF11; plasmolipin
Accession No.	Swiss-Prot#:Q9Y342NCBI Gene ID:51090
Uniprot	Q9Y342
GeneID	51090;
Calculated MW	20kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C


Application Details

WB 1:500 - 1:2000IF 1:50 - 1:200

Images



Western blot analysis of extracts of various cell lines, using PLL_P at 1:1000 dilution.



Immunofluorescence analysis of L929 cells using PLLP at dilution of 1:100. Blue: DAPI for nuclear staining.

Background

Note: This product is for in vitro research use only