

MAP7D1 Polyclonal Antibody

Catalog No: #29458



Package Size: #29458-1 50ul #29458-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	MAP7D1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Mouse
Immunogen Description	Recombinant fusion protein of human MAP7D1 (NP_060537.3).
Other Names	MAP7D1; PARCC1; RPRC1; MAP7 domain containing 1
Accession No.	Swiss-Prot#:Q3KQU3NCBI Gene ID:55700
Uniprot	Q3KQU3
GeneID	55700;
Calculated MW	110kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using MAP7D1 at 1:1000 dilution.

Background

Note: This product is for in vitro research use only