

# ITFG1 Polyclonal Antibody

Catalog No: #29482



Package Size: #29482-1 50ul #29482-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

|                       |                                                                           |
|-----------------------|---------------------------------------------------------------------------|
| Product Name          | ITFG1 Polyclonal Antibody                                                 |
| Host Species          | Rabbit                                                                    |
| Clonality             | Polyclonal                                                                |
| Isotype               | IgG                                                                       |
| Purification          | Affinity purification                                                     |
| Applications          | WB                                                                        |
| Species Reactivity    | Human,Mouse,Rat                                                           |
| Immunogen Description | Recombinant fusion protein of human ITFG1 (NP_110417.2).                  |
| Other Names           | ITFG1; 2310047C21Rik; CDA08; LNKN-1; TIP; T-cell immunomodulatory protein |
| Accession No.         | Swiss-Prot#:Q8TB96NCBI Gene ID:81533                                      |
| Uniprot               | Q8TB96                                                                    |
| GeneID                | 81533;                                                                    |
| Calculated MW         | 80kDa                                                                     |
| Formulation           | Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.         |
| Storage               | Store at -20°C                                                            |

## Application Details

WB 1:500 - 1:2000

## Images



Western blot analysis of extracts of various cell lines, using ITFG1 at 1:1000 dilution.

## Background

---

Note: This product is for in vitro research use only