

CD99L2 Polyclonal Antibody

Catalog No: #29486

Package Size: #29486-1 50ul #29486-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	CD99L2 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of human CD99L2 (NP_113650.2).
Other Names	CD99L2; CD99B; MIC2L1; CD99 molecule like 2
Accession No.	Swiss-Prot#:Q8TCZ2NCBI Gene ID:83692
Uniprot	Q8TCZ2
GeneID	83692;
Calculated MW	28kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using CD99L2 at 1:1000 dilution.

Background

This gene encodes a cell-surface protein that is similar to CD99. A similar protein in mouse functions as an adhesion molecule during leukocyte extravasation. Alternate splicing results in multiple transcript variants.

Note: This product is for in vitro research use only