

ZNF766 Polyclonal Antibody

Catalog No: #29494

Package Size: #29494-1 50ul #29494-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	ZNF766 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Mouse,Rat
Immunogen Description	Recombinant fusion protein of human ZNF766 (NP_001010851.1).
Other Names	ZNF766; zinc finger protein 766
Accession No.	Swiss-Prot#:Q5HY98NCBI Gene ID:90321
Uniprot	Q5HY98
GeneID	90321;
Calculated MW	60kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using ZNF766 at 1000 dilution.

Background

Note: This product is for in vitro research use only