

MVB12A Polyclonal Antibody

Catalog No: #29497



Package Size: #29497-1 50ul #29497-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	MVB12A Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Mouse
Immunogen Description	Recombinant fusion protein of human MVB12A (NP_612410.1).
Other Names	MVB12A; CFBP; FAM125A; multivesicular body subunit 12A
Accession No.	Swiss-Prot#:Q96EY5NCBI Gene ID:93343
Uniprot	Q96EY5
GeneID	93343;
Calculated MW	29kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of mouse liver, using MVB12A at 1000 dilution.

Background

Note: This product is for in vitro research use only