

TMEM44 Polyclonal Antibody

Catalog No: #29498

Package Size: #29498-1 50ul #29498-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	TMEM44 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human
Immunogen Description	Recombinant fusion protein of human TMEM44 (NP_001159777.1).
Other Names	TMEM44; transmembrane protein 44
Accession No.	Swiss-Prot#:Q2T9K0NCBI Gene ID:93109
Uniprot	Q2T9K0
GeneID	93109;
Calculated MW	52kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of A-431 cells, using TMEM44 at 1:1000 dilution.

Background

Note: This product is for in vitro research use only