

RBP7 Polyclonal Antibody

Catalog No: #29504

Package Size: #29504-1 50ul #29504-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	RBP7 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Mouse,Rat
Immunogen Description	Recombinant fusion protein of human RBP7 (NP_443192.1).
Other Names	RBP7; CRABP4; CRBP4; CRBP4V; retinol binding protein 7
Accession No.	Swiss-Prot#:Q96R05NCBI Gene ID:116362
Uniprot	Q96R05
GeneID	116362;
Calculated MW	15kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using RBP7 at 1000 dilution.

Background

The protein encoded by this gene is a member of the cellular retinol-binding protein (CRBP) family, whose members are required for vitamin A stability and metabolism. The encoded protein binds all-trans-retinol and is structurally similar to other CRBPs; however, it has a lower binding affinity for retinol than other CRBPs.

Note: This product is for in vitro research use only