

CHCHD1 Polyclonal Antibody

Catalog No: #29506



Package Size: #29506-1 50ul #29506-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	CHCHD1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Mouse
Immunogen Description	Recombinant fusion protein of human CHCHD1 (NP_976043.1).
Other Names	CHCHD1; C10orf34; C2360; MRP-S37; coiled-coil-helix-coiled-coil-helix domain containing 1
Accession No.	Swiss-Prot#:Q96BP2NCBI Gene ID:118487
Uniprot	Q96BP2
GeneID	118487;
Calculated MW	13kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using CHCHD1 at 1000 dilution.

Background

Note: This product is for in vitro research use only