

TARSL2 Polyclonal Antibody

Catalog No: #29507



Package Size: #29507-1 50ul #29507-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	TARSL2 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of human TARSL2 (NP_689547.2).
Other Names	TARSL2; threonyl-tRNA synthetase like 2
Accession No.	Swiss-Prot#:A2RTX5NCBI Gene ID:123283
Uniprot	A2RTX5
GeneID	123283;
Calculated MW	110kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using TARSL2 at 1:1000 dilution.

Background

Note: This product is for in vitro research use only