

LSM12 Polyclonal Antibody

Catalog No: #29508



Package Size: #29508-1 50ul #29508-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	LSM12 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Mouse
Immunogen Description	Recombinant fusion protein of human LSM12 (NP_689557.1).
Other Names	LSM12; PNAS-135; LSM12 homolog
Accession No.	Swiss-Prot#:Q3MHD2NCBI Gene ID:124801
Uniprot	Q3MHD2
GeneID	124801;
Calculated MW	22kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using LSM12 at 1000 dilution.

Background

Note: This product is for in vitro research use only