

TPRG1L Polyclonal Antibody

Catalog No: #29510



Package Size: #29510-1 50ul #29510-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	TPRG1L Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Mouse,Rat
Immunogen Description	Recombinant fusion protein of human TPRG1L (NP_877429.2).
Other Names	TPRG1L; FAM79A; h-mover; mover; tumor protein p63 regulated 1 like
Accession No.	Swiss-Prot#:Q5T0D9NCBI Gene ID:127262
Uniprot	Q5T0D9
GeneID	127262;
Calculated MW	30kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using TPRG1L at 1:1000 dilution.

Background

Note: This product is for in vitro research use only