

ACTL9 Polyclonal Antibody

Catalog No: #29523



Package Size: #29523-1 50ul #29523-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	ACTL9 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Mouse,Rat
Immunogen Description	Recombinant fusion protein of human ACTL9 (NP_848620.3).
Other Names	ACTL9; HSD21; actin like 9
Accession No.	Swiss-Prot#:Q8TC94NCBI Gene ID:284382
Uniprot	Q8TC94
GeneID	284382;
Calculated MW	46kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using ACTL9 at 1:1000 dilution.

Background

Note: This product is for in vitro research use only