

RFPL4 Polyclonal Antibody

Catalog No: #29526



Package Size: #29526-1 50ul #29526-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	RFPL4 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Mouse,Rat
Immunogen Description	A synthetic peptide of human RFPL4 (NP_001138486.1).
Other Names	RFPL4A; RFPL4; RNF210; ret finger protein-like 4A
Accession No.	Swiss-Prot#:A6NLU0NCBI Gene ID:342931
Uniprot	A6NLU0
GeneID	342931;
Calculated MW	32kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using RFPL4 at 1:1000 dilution.

Background

Note: This product is for in vitro research use only