

## N4-Methylcytidine Polyclonal Antibody

Catalog No: #29530

Package Size: #29530-1 50ul #29530-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	N4-Methylcytidine Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IF
Species Reactivity	Human
Immunogen Description	A synthetic peptide of human 4mC
Accession No.	Swiss-Prot#:_NCBI Gene ID:_
Calculated MW	Refer to figures
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

## Application Details

WB□1:500 - 1:2000IF□1:50 - 1:200

## Images



Immunofluorescence analysis of HeLa cells using 4mC at dilution of 1:100 (40x lens). Blue: DAPI for nuclear staining.

## Background

Note: This product is for in vitro research use only