

ZNRF3 Polyclonal Antibody

Catalog No: #29540



Package Size: #29540-1 50ul #29540-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	ZNRF3 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IHC
Species Reactivity	Human,Mouse,Rat
Immunogen Description	A synthetic peptide of human ZNRF3 (NP_115549.2).
Other Names	ZNRF3; BK747E2.3; RNF203; zinc and ring finger 3
Accession No.	Swiss-Prot#:Q9ULT6NCBI Gene ID:84133
Uniprot	Q9ULT6
GeneID	84133;
Calculated MW	120kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C


Application Details

WB 1:500 - 1:2000IHC 1:50 - 1:200


Images




Western blot analysis of extracts of various cell lines, using ZNRF3 at 1:1000 dilution.




Immunohistochemistry of paraffin-embedded human colon carcinoma using ZNRF3 at dilution of 1:100 (40x lens).




Immunohistochemistry of paraffin-embedded human colon using ZNRF3 at dilution of 1:100 (40x lens).



Immunohistochemistry of paraffin-embedded human liver using ZNRF3 at dilution of 1:100 (40x lens).



Immunohistochemistry of paraffin-embedded mouse spleen using ZNRF3 at dilution of 1:100 (40x lens).



Immunohistochemistry of paraffin-embedded rat lung using ZNRF3 at dilution of 1:100 (40x lens).

Note: This product is for in vitro research use only