

C1QTNF7 Polyclonal Antibody

Catalog No: #29543



Package Size: #29543-1 50ul #29543-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	C1QTNF7 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Mouse,Rat
Immunogen Description	Recombinant fusion protein of human C1QTNF7 (NP_001128642.1).
Other Names	C1QTNF7; CTRP7; ZACRP7; C1q and TNF related 7
Accession No.	Swiss-Prot#:Q9BXJ2NCBI Gene ID:114905
Uniprot	Q9BXJ2
GeneID	114905;
Calculated MW	Refer to figures
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB□1:500 - 1:2000

Background

Note: This product is for in vitro research use only